



Evidence Series:

**Discrepancy between
prescribed and actual APD
prescription delivery:**
Identification using cyclor
remote management
technology

Catherine A Firanek et al

Firanek C et al, MP557 Discrepancy between prescribed and actual APD prescription delivery: Identification using cyclor remote management technology Neph Dial Trans 2017; 32 (suppl 3): iii633

Vantive, Homechoice Claria, and Sharesource are trademarks of Vantive Health LLC or its affiliates.

Sharesource
REMOTE PATIENT MONITORING



BACKGROUND

Historically, clinicians have been unable to proactively identify patients missing or shortening PD treatments.

NON-ADHERENCE TO
>10% of the Peritoneal
Dialysis (PD)
prescription is
associated with

TECHNIQUE FAILURE
PERITONITIS
HOSPITALISATIONS
& MORTALITY^{1,2}

1. J Bernardini, M Nagy, B Piraino. Pattern of Noncompliance with Dialysis Exchanges in Peritoneal Dialysis Patients. Am J Kidney Dis 2000; 35: 1104-1110.

2. J Bernardini, B Piraino. Compliance in CAPD and CCPD Patients as Measured by Supply Inventories During Home Visits. Am J Kidney Dis 1998; 31: 107-107.



**Automated
Peritoneal
Dialysis
(APD)**

cyclers embedded with Remote Patient Management (RPM) technology can detect early treatment-related issues, allowing intervention to potentially prevent clinically significant events.



OBJECTIVES

- > Evaluate actual APD treatment time compared with prescribed treatment time, using an APD device with embedded RPM technology (Homechoice Claria with Sharesource).
- > Determine if Clinicians using APD with Sharesource have greater visibility to patient adherence patterns to allow early intervention.



ENDPOINTS

Patient adherence and early intervention



METHODS

DATA ON
399 European APD
PATIENTS
WERE ANALYSED

- > Patients with ≥ 3 months on the Homechoice Claria with Sharesource were examined for weekly treatment frequency and actual versus prescribed treatment time.
- > An assumption was made that patients perform APD therapy 7 days per week.
- > Patients with gaps in treatment >30 days were omitted/excluded.

- > Any treatments occurring in the first 14 days from the very first available treatment were considered as training time and were excluded.
- > Time (days) on treatment was determined from the first treatment after the training period to the last available treatment for a patient.
- > Weekly treatment frequency was the number of treatments in Sharesource/30 x 7. Eg., if a patient had 27 treatments out of 30 days, then Weekly rate = $(27/30) \times 7 = 6.3$.
- > Treatment differences were treatment time prescribed – actual treatment time performed.

RESULTS

During the 1st month of therapy:


- > 30% (115) of patients missed ≥ 4 treatments (>10% of prescribed therapy)
- > 12% (47) of patients missed ≥ 9 treatments

In the 1st week of therapy:


- > 24.3% (97) and 9.5% (38) of patients had ≥ 100 minutes and ≥ 300 minutes, respectively, less actual therapy time than prescribed



Number of Missed Treatments in First Month of Dialysis (N=399)



Number of Patients Who Missed Significant Treatment Time / Week by Week of Therapy (N=399)





RESULTS

In combined results of weeks 1, 2, 4 and 6:

> 43% of patients missed >5% of prescribed dwell time


20.6% OF PATIENTS
MISSED >10%

OF PRESCRIBED DWELL TIME

> 11.9% missed >15% of prescribed dwell time

> 7% missed >20% of prescribed dwell time

Number of Patients With at Least 3 Months of Treatment Who Missed Significant Dwell Time by Week of Therapy (N=399)



CONCLUSIONS

- > Current standard of care does not allow visibility to determine adherence to prescribed PD therapy.
- > Sharesource remote patient management platform allows clinicians to securely view their patients' daily home dialysis treatment data.
- > Visibility to adherence patterns may provide opportunities for clinicians to intervene, educate or retrain the patient in a more timely manner.

- > Clinicians using APD with Sharesource have greater visibility to

**PATIENT ADHERENCE
PATTERNS
WHICH MAY ALLOW
EARLY
INTERVENTION**

